7/28

- Optimize the condition of folding DNA origami – the mole ratio between scaffold and staple.

Mole ratio(scaffold:staple)	scaffold only	staple only	1:5	1:7	1:10
scafold(100nM)	52	0	2	2	2
staple(250nM)	0	8	4	5.6	8
10X TE	2	2	2	2	2
MgCl ₂ (200mM)	2	2	2	2	2
ddH ₂ O	14	8	10	8.4	6
total	20	820	20	20	20

7/29

- Run the gel electrophoresis
 - Run two gels for each of the three structure. Six gels in total.
 - Box: 1-264(**A**)
 - bottom:265-348(**B**)
 - Top : 349-434(**C**)
 - 7 wells on one gel with the first one and the last one being the ladder.
Load 20 micro l in each well.
 - Choose the best ratio based on the gel result and conduct the PCR under different Mg²⁺ ion concentration accordingly.

Side notes:


- PCR notation: f for scaffold, p for staple, Box: 1-264 are noted as **A**, bottom:265-348 are noted as (**B**)
- scaffold, staple (note 1, 256, 349 respectively),200mMMgCl₂, are all put in the -20 degree celcius fridge.

- Steps for running gel electrophoresis**


1. Mix 1xTAE buffer 150mL and agarose gel powder 1.5g, and heat up the mixture with microwave till the powder is completely dissolved and the solution becomes transparent.
2. Wait till the solution cools down a bit. Pour the mixture into the glue plate and put on the suitable tooth comb; leave it for 50 minutes for solidification.
3. Put the solid agar gel with the plastic case into the electrophoresis groove.
4. Mix 20 μ L of sample solution and 2.2 μ L of 10X loading dye, and load each sample to wells. (the first well and the last well is loaded with DNA markers.)
5. Add 3 μ L of (10^-5M) EtBr solution to the electrophoresis groove evenly.
6. Put on the lid and open the switch to run the electrophoresis with 110V half-voltage for an hour. (cover with tin foil to prevent EtBr from decomposing)
7. After the electrophoresis, put the agar gel in the UV box for observation.

- **determining scaffold/staple ratio**


- box result (marker/scaffold/staple/1:5/1:5/1:7/1:7/1:10/1:10/marker)
 - We can see that all lanes have folded products, but we are not sure which position is the correct product. It seems that 1:5 and 1:10 are the two most productive ratios.

- bottom lid result
(marker/scaffold/staple/1:5/1:5/1:7/1:7/1:10/1:10/marker)

- top lid result (marker/scaffold/staple/1:5/1:5/1:7/1:7/1:10/1:10/marker)

7/31

- Test 5 different Mg²⁺ ion concentrations. (10,15,20,25,30 microM)
 - $5 \times 2 + 1 + 1 + 2 = 14$ wells in total. 3 gels in total.
 - Since the heat released from the gel electrophoresis will disentangle the structure, we put the TAE buffer used during the experiment under 4°C for 30 min.

f&p formula. 20 micro l in total for each.

The image shows a handwritten lab notebook page. At the top, it says "30 min 10% ac". Below that is a sketch of a gel electrophoresis tank with two lanes labeled "f" and "p". To the left of the gel, there is a table for preparing "f" and "p" formulas. The table has two columns: "f" and "p". Under "f", there are four rows with values: 2, 2, 0, and 0. Under "p", there are four rows with values: 2, 3, 0, and 4. To the left of the table, there are handwritten notes: "10X TE 2", "MgCl₂ 30", "f 2", "p 4", and "H₂O 9".


	f	p
	2	2
	3	3
	0	0
	0	4
H ₂ O	9	11

remember to run PCR for f & p as well.

Side notes :

1. The red dots on the staple are used.
2. Marker needs to be put in the -20 degree celsius firdge.
3. gel electrophoresis takes about 30 min. Wait till the deep blue dye reaches the third to last line.


Result :

Box - choose 10mM
(marker/scaffold/staple/10mM/10mM/15mM/15mM/20mM/20mM/25mM/25mM/30mM
/30mM/marker)


8/1

- Run the gel once again. The same as yesterday.

bottom lid - choose 15mM

(marker/scaffold/staple/10mM/10mM/15mM/15mM/20mM/20mM/25mM/25mM/30mM
/30mM/marker)

top lid - choose 15mM


(marker/scaffold/staple/10mM/10mM/15mM/15mM/20mM/20mM/25mM/25mM/30mM /30mM/marker)

- box(10mM MgCl₂), bottom top(15mM MgCl₂) 50 μl for 10 tubes in total.
Run the PCR **using NanoMuscle protocol.**

8/2

- **DNA origami Purification**
***wash before purify**
 1. assemble the filter and centrifuge tube, load 200μl 1xTEbuffer with 10mM Mg²⁺(box).
 2. 200μl 1xTEbuffer with 15mM Mg²⁺(lid)
 3. spin down at 14000g for 3 min
 - a. **7000g is suggested since little unbalance could lead to more serious result if the rotation is as quick as 14000g.**
 4. Turn the filter upside. (no need to cap). Spin the device at **1000 g for 2 min.**

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※wash filter

1. wash with ddH2O*2
 - a. rotate in the 'up position' at 7000g for 3min.
 - b. rotate upside at 100g for 2 min.
2. sonication in ddH2O for 30~60min
3. keep in 20%EtOH

※plan for 8/4


1. run three gels, add 2micro l loading dye of 6x conc. 10 micro l for each well.
2. original sample/U1/D1/U2/D2/U3/D3/U4/D4/U5/D5(x and t look quite similar, please pay extra attention to them.)

8/4

- Run the gel after purification

bottom lid (B) :


[from left to right] ladder/ref./U1/U2/U3/U4/U5/D1/D2/D3/D4/D5/ladder

- D is a bit weird. Nothing is seen.
- Purification is pretty good (U is apparently less in the third round). Maybe purified for three times is proper.

top lid (t):

[from left to right] ladder/ref./U1/U2/U3/U4/U5/D1/D2/D3/D4/D5/ladder

- Purification is pretty good. Purify 4 times is proper.

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1. Pick up a TEM copper grid (small) with negative tweezer.
2. Transfer 7 uL sample solution onto the TEM copper grid.
3. Wait for 35 seconds.
4. Remove sample solution from the side of the droplet with the edge of the filter paper.
5. Transfer a droplet of 2% UA (negative stain solution) onto the TEM copper grid.
6. Wait for 30 seconds.
7. Remove negative stain solution from the side of the droplet with the edge of the filter paper.
8. Place TEM copper grid under a lamp for 20 minutes.
9. Put TEM copper grid into TEM. (assisted by technician)
10. Observe with TEM and save the coordinates you intend to take photo of.
11. Take pictures and save files. (assisted by technician)

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
10/1

1. Add blocking staples solution to base solution. The mole of blocking staples is twenty times as the scaffold.
2. Put the mixture of base and blocker in Mini-Shaking-Hybridisation Oven, and shake at 25 °C for one hour.
1. Add hinge staples solution to lid solution. The mole of hinge staples is five times as the scaffold.
2. Put the mixture of lid and hinge in Mini-Shaking-Hybridisation Oven, and shake at 25 °C for one hour.
3. Purify the solution of lid with hinge to remove excess staples.

10/2

1. Mix the blocked base and purified left lid with hinge and purified right lid with hinge
2. Put the mixture in Mini-Shaking-Hybridisation Oven, and shake at 30 °C overnight.
3. Purify the solution of assemblies to remove excess staples.

10/3



sample assembly and purification

Corresponding to the previous results, both Right and Left Lids exist between 4kb and 5kb, and the base appeared near 5kb. Notably, the components weren't purified in these three lanes. Since blockers and hinges are all short DNA fragments, we expect that no obvious band shifts will happen upon adding them, which matches the experiment results. In Figure 10 we could also observe that for Left and Right Lids, the bands of excess DNA slightly faded more than the band of Left/Right Lids + Hinges after the purifying process, so we decided to purify both Right/Left Lid + Hinge before assembling them with the base. Finally, the products are barely seen in the lane of Assembled**, so we decided to purify the Assembled components once only.

10/4

sonicate for ten minutes

BTX

10/5

1. To flip left lid: Add signal L to the mixture of Base and Both Lid
- To flip right lid: Add signal R to the mixture of Base and Both Lid
- To flip Both lids: Add signal L and signal R to the mixture of Base and Both Lid
(add signal according to the fold below)

signal:blocker	BTX	signalL	signalR
10:1	60	12	
	60		12
	60	12	12
7:1	30	4.2	
	30		4.2
	30	4.2	4.2
5:1	30	3	
	30		3
	30	3	3


2. Put the mixture in Mini-Shaking-Hybridisation Oven, and shake at 30 °C overnight.

10/8


flip the left lid(BTX-FL)

signal:blocker=10:1

sonicate for 10min

the shape is similar with FL, but the size is two time bigger than our structure


the shape and size is similar to flip the left lid(BTX-FL)

10/9

flip the right lid(BTX-FR)

signal:blocker=10:1

sonicate for 10min

10/14

flip both lids

the lids didn't attach to the base properly